


16 Priority PAHs and Methylnaphthalenes on Rxi-5Sil MS using the GC Accelerator Kit in a 120 V Oven



| Peaks | tR (min) | Conc. (µg/mL) |
|-------------------------------|----------|---------------|
| 1. Naphthalene-d8 (IS) | 3.44 | 2.0 |
| 2. Naphthalene | 3.45 | 1.0 |
| 3. 2-Methylnaphthalene | 3.86 | 1.0 |
| 4. 1-Methylnaphthalene | 3.92 | 1.0 |
| 5. 2-Fluorobiphenyl (SS) | 4.08 | 2.0 |
| 6. Acenaphthylene | 4.41 | 1.0 |
| 7. Acenaphthene-d10 (IS) | 4.49 | 2.0 |
| 8. Acenaphthene | 4.51 | 1.0 |
| 9. Fluorene | 4.83 | 1.0 |
| 10. 2,4,6-Tribromophenol (SS) | 4.97 | 2.0 |
| 11. Phenanthrene-d10 (IS) | 5.40 | 2.0 |
| 12. Phenanthrene | 5.41 | 1.0 |
| 13. Anthracene | 5.45 | 1.0 |
| 14. Fluoranthene | 6.15 | 1.0 |
| 15. Pyrene-d10 (SS) | 6.28 | 2.0 |
| 16. Pyrene | 6.29 | 1.0 |
| 17. p-Terphenyl-D14 (SS) | 6.37 | 2.0 |
| 18. Benz[a]anthracene | 7.12 | 1.0 |
| 19. Chrysene-d12 (IS) | 7.13 | 2.0 |
| 20. Chrysene | 7.15 | 1.0 |
| 21. Benzo[b]fluoranthene | 8.17 | 1.0 |
| 22. Benzo[k]fluoranthene | 8.21 | 1.0 |
| 23. Benzo[a]pyrene | 8.56 | 1.0 |
| 24. Perylene-d12 (IS) | 8.64 | 2.0 |
| 25. Indeno[1,2,3-cd]pyrene | 10.16 | 1.0 |
| 26. Dibenzo[a,h]anthracene | 10.20 | 1.0 |
| 27. Benzo[ghi]perylene | 10.58 | 1.0 |

Column Rxi-5Sil MS, 20 m, 0.15 mm ID, 0.15 µm (cat.# 43816)
Sample EPA Method 8310 PAH mixture (cat.# 31874)
 Revised SV internal standard mix (cat.# 31886)
 Revised B/N surrogate mix (cat.# 31888)
 Acid surrogate mix (4/89 SOW) (cat.# 31063)
 Dichloromethane

Diluent:
Injection
 Inj. Vol.: 1.0 µL split (split ratio 20:1)
 Liner: Topaz 4 mm single taper w/wool (cat.# 23303)
 Inj. Temp.: 275 °C

Oven
 Oven Temp.: 60 °C (hold 0.7 min) to 285 °C at 39.8 °C/min to 305 °C at 4.3 °C/min to 320 °C at 28.5 °C/min (hold 3.5 min)

Carrier Gas He, constant flow
Flow Rate: 1.0 mL/min
Detector MS
Mode: SIM

| Group | Start Time (min) | Ion(s) (m/z) | Dwell (ms) |
|-------|------------------|--|------------|
| 1 | 3.29 | 102, 108, 128, 136 | 25 |
| 2 | 3.71 | 85, 115, 142.1, 172.1 | 20 |
| 3 | 4.28 | 76, 82, 152.1, 153.1, 164.1 | 20 |
| 4 | 4.71 | 82.40, 142.90, 166.1, 329.8 | 25 |
| 5 | 5.24 | 89, 94, 178.1, 188.1 | 15 |
| 6 | 5.88 | 101, 106.1, 122.1, 202.1, 212.1, 244.1 | 20 |
| 7 | 6.83 | 101, 120.1, 226.1, 228.1, 240.2 | 10 |
| 8 | 7.77 | 126, 252.1 | 25 |
| 9 | 8.44 | 126, 132.1, 252.1, 264.2 | 25 |
| 10 | 9.57 | 138, 139, 276.1, 278.1 | 25 |

Transfer Line Temp.: 280 °C
 Analyzer Type: Quadrupole
 Source Type: Extractor
 Extractor Lens: 9 mm ID
 Source Temp.: 330 °C
 Quad Temp.: 180 °C
 Solvent Delay Time: 1 min
 Tune Type: DFTPP
 Ionization Mode: EI
Instrument Agilent 7890B GC & 5977A MSD
Notes Fast SIM analysis of 16 priority PAHs plus the methylnaphthalenes in a 120 V oven equipped with the GC Accelerator kit (cat.# 23849) (injected 1 µg/mL = 0.05 ng on-column).